

TECHNICAL SPECIFICATIONS

BHQ-1 Quencher - Chip

Product information

	BHQ-1 Quencher - Chip
Properties	Single-use; microcolumn incl.
Chip Type	Microfluidic
End Product	5'-6-FAM / 3'-BHQ-1 labeled DNA oligo
Material	Polymer
Shelf Life	6 months
Short-term Storage	Dry; 2-8°C; dark
Long-term storage	Dry; 2-8°C; dark
Combination	6-FAM DNA - Cartridge
Scope of Delivery	1 microfluidic chip; OliPure HIC Purification Kit
Required	6-FAM DNA - Cartridge; Kilobaser one-XT