

TECHNICAL SPECIFICATIONS

Phosphorylation - Chip

Product information

	Phosphorylation - Chip
Properties	Single-use; microcolumn incl.
Chip Type	Microfluidic
End Product	5'-phosphorylated DNA oligo
Material	Polymer
Shelf Life	1 year
Short-term Storage	Dry; room temperature
Long-term storage	Dry; room temperature
Combination	Phosphorylation DNA - Cartridge
Scope of Delivery	1 microfluidic chip
Required	Phosphorylation DNA - Cartridge; Kilobaser one-XT