

# TECHNICAL SPECIFICATIONS

## Standard - Chip

### Product information

	<b>Standard - Chip</b>
<b>Properties</b>	Single-use; microcolumn incl.
<b>Chip Type</b>	Microfluidic
<b>End Product</b>	Unlabeled DNA oligo
<b>Material</b>	Polymer
<b>Shelf Life</b>	1 year
<b>Short-term Storage</b>	Dry; room temperature
<b>Long-term storage</b>	Dry; room temperature
<b>Combination</b>	All Kilobaser cartridges
<b>Scope of Delivery</b>	1 microfluidic chip
<b>Required</b>	Any cartridge; Kilobaser one or Kilobaser one-XT